



Running the Quant-iT assay on the iD3 plate reader

Adapted from manufacturer's protocol.

Updated 29-Feb-2024 by Katie Murphy

HELPFUL HINTS:

- If your kit includes separate dye and buffer, allow the kit to come to room temperature before use (the dye needs to be thawed). This will take at least 30 minutes.
- Each assay plate needs to include the 8 standards. This means you can run up to 88 of your samples per plate.
- You do NOT need to run full plates. Leave any wells empty that you don't need (don't waste reagents please).
- Use solid black plates with flat-bottom wells.
- Turn on the plate reader prior to setting up your assay plate, so it'll be ready to go when needed.
- For more accurate pipetting and therefore more accurate quantification, we recommend using at least 2 μ L DNA rather than just 1 μ L.
- For your calculations, it is simplest if you use the same volume of DNA for all samples on a plate. But this isn't strictly necessary, as long as you record what volume of DNA you used for each sample.

MATERIALS & REAGENTS:

- DNA to be quantified
- Reservoir (for adding assay buffer to plate with a multichannel pipette)
- HS or BR Quant-iT kit
 - HS is recommended if you expect samples to be < 100 ng/ μ L
- Solid black 96-well plate (Corning #3915)
- Foil plate seal



PROTOCOL:

Step 1 – Setting up your assay plate

Please follow 1A if you plan to use the pre-defined software template (2A). Alternatively, you can follow 1B if you aren't using the pre-defined software protocol, which gives you flexibility in where you run your standards.

1A. Setting up your assay plate to use with our pre-defined software template

If using this protocol, you **MUST** load your standards in the 12th column of your assay plate.

- 1A.1 If not using a 1X kit: you need to make a 1X working stock. Calculate the number of reactions needed, making sure to include a few extra reactions to allow for pipetting error. For each reaction, you'll need 1 μ L of dye and 200 μ L of buffer.
 - a. Example: for a full plate of 96 reactions, mix enough for 100 reactions – you'd need 100 μ L of dye and 20mL of buffer.
 - b. Once you dilute the dye, it's stable at room temp for 3 hours as long as it's protected from light.
- 1A.2 Add 200 μ L of diluted dye to each well of your 96 well plate that you plan to use for samples or standards.
- 1A.3 Lightly mix your standards by inverting or vortexing, then briefly spin down. Add 10 μ L of each standard to the following wells in column 12.
 - a. A12 – standard 1 (least concentrated)
 - b. B12 – standard 2
 - c. C12 – standard 3
 - d. D12 – standard 4
 - e. E12 – standard 5
 - f. F12 – standard 6
 - g. G12 – standard 7
 - h. H12 – standard 8 (most concentrated)
- 1A.4 Add 1-20 μ L of each DNA sample to its designated well of the assay plate.
 - a. We recommend 2 μ L instead of 1 μ L if possible, as pipetting is more accurate with 2 μ L.
- 1A.5 Seal the assay plate thoroughly with a foil seal, then vortex to mix.
- 1A.6 Spin plate briefly to collect contents to the bottom of the wells.
- 1A.7 Incubate plate at room temperature for at least 2 minutes before reading on the iD3 plate reader.
 - a. Your assay will be stable for up to two hours at room temperature, as long as you store the plate in a dark location (or keep covered with foil).
- 1A.8 Remove the foil seal before inserting your plate into the plate reader.

1B. Setting up your assay plate without using our pre-defined software protocol.

If using this protocol, you can load your standards anywhere on your assay plate – even in a separate plate. If you have a full plate of libraries, this would allow you to run all samples in one plate, and the standards in a second plate. Please note, you must follow Step 2B if using this protocol.

- 1B.1 If not using a 1X kit: you need to make a 1X working stock. Calculate the number of reactions needed, making sure to include a few extra reactions to allow for pipetting error. For each reaction, you'll need 1 μ L of dye and 200 μ L of buffer.
 - a. Example: for a full plate of 96 reactions, mix enough for 100 reactions – you'd need 100 μ L of dye and 20mL of buffer.
 - b. Once you dilute the dye, it's stable at room temp for 3 hours as long as it's protected from light.



- 1B.2 Add 200 μ L of diluted dye to each well of your 96 well plate(s) that you plan to use for samples or standards.
- 1B.3 Lightly mix your standards by inverting or vortexing, then briefly spin down. Add 10 μ L of each standard to the wells of your choice. Make sure to record which wells you use.
- 1B.4 Add 1-20 μ L of each DNA sample to its designated well of the assay plate(s).
 - a. We recommend 2 μ L instead of 1 μ L if possible, as pipetting is more accurate with 2 μ L.
- 1B.5 Seal the assay plate(s) thoroughly with a foil seal, then vortex to mix.
- 1B.6 Spin plate(s) briefly to collect contents to the bottom of the wells.
- 1B.7 Incubate plate(s) at room temperature for at least 2 minutes before reading on the iD3 plate reader.
 - a. Your assay will be stable for up to two hours at room temperature, as long as you store the plate in a dark location (or keep covered with foil).
- 1B.8 Remove the foil seal before inserting your plate into the plate reader.

Step 2 – Running your quantification assay using the iD3 plate reader

You can either use a pre-defined protocol (2A), or you can set up your own run from scratch (2B).

2A. Using LAB's pre-defined plate reader protocol & analyzing within SoftMax Pro

If using this protocol, you MUST have set up your standards in column 12 of each assay plate.

- 2A.1 If not already on, turn on both the plate reader and the computer.
- 2A.2 Log onto the computer using our standard lab login info (it's printed on the side of the computer if you forget).
- 2A.3 Click on the corresponding desktop icon for your assay.
 - a. If you're running the High Sensitivity kit, click on "HS Quant-iT assay" icon.
 - b. If you're using the Broad Range kit, click on "BR Quant-iT kit" icon.
 - c. *Make sure you open the correct protocol for your assay type. If you open the wrong one, your standard curve will be incorrect, and therefore your concentrations will be incorrect.*
- 2A.4 Verify that the software is connected to the plate reader. On the upper left of the window, you should see a green checkmark on the SpectraMax device icon. If instead there is a red circle with a slash through it, do the following:
 - a. Double-click the SpectraMax icon.
 - b. After the window pops up and populates with a list of instruments, click on the SpectraMax iD3 and connect to it. It should be the only instrument listed.
 - c. Close this pop-up window if it doesn't automatically do so.
- 2A.5 Open the plate reader drawer by clicking on the "Open/Close" button in the menu bar of the software, or by clicking the eject button on the iD3 touchscreen.
- 2A.6 Insert your plate, with A1 on the upper left and H12 on the bottom right.
- 2A.7 Close the drawer by clicking on the "Open/Close" button in the menu bar of the software, or by clicking the eject button on the iD3 touchscreen.
 - a. *DO NOT MANUALLY PUSH IN THE DRAWER BY HAND!!!*
- 2A.8 Click the "Read" button (green arrow).
- 2A.9 Click "skip/okay" if there is a popup about optimization (you don't need to check any boxes).
- 2A.10 The read will take approximately 1 min 30 seconds per plate.
- 2A.11 Once the read is complete, verify that your standard curve looks correct. Remove/mask any incorrect standard value(s).
- 2A.12 We highly recommend forcing the y-intercept to be the fluorescent value obtained for your first standard (0 ng). You can do this in the SoftMax Pro software by following these steps:
 - a. Scroll until you are viewing the standard curve.



- b. From the buttons that are now along the top of the software screen, select “Curve Fit Settings”.
 - c. In the pop-up that appears, set value “A (the y-intercept) to the fluorescent value obtained for standard 1.
 - d. Save/close the pop-up menu.
 - e. Your standard curve should now show an adjusted linear equation.
 - f. You can verify you’ve done this correctly by scrolling to the table with your standard values. The back-calculated value for the first standard should now be exactly 0 ng.
- 2A.13 Record the total ng (“Result” column) reported for your samples.
- a. You can copy/paste this column into a spreadsheet, or you can simply export the entire output to excel.
 - b. To export your data in Excel (.xls) format:
 - i. Click on the Plate icon in the top left of the screen.
 - ii. Select Export.
 - iii. Select your entire experiment (or plates of interest) & export your raw data.
 - iv. Select “column” or “plate” format, depending on your preference.
 - v. Change file type to .xls.
- 2A.14 Divide your result by the volume (μL) of DNA you used for each sample to get your original ng/ μL concentration.
- a. For example, if you only used 1 μL of DNA for your samples, then the outputted result is your ng/ μL .
 - b. If you used 5 μL of DNA, then divide the values by 5 to get your original ng/ μL concentration.

2B. Setting up the plate reader from scratch & analyzing in Excel.

Use this protocol if you set up your standards in wells other than column 12, or if you prefer to analyze your data in excel yourself.

- 2B.1 If not already on, turn on both the plate reader and the computer.
- 2B.2 Log on to the computer using our standard lab login info (it’s printed on the side of the computer if you forget).
- 2B.3 Launch the SoftMax Pro software by double-clicking the desktop icon.
- 2B.4 Verify that the software is connected to the plate reader. On the upper left of the window, you should see a green checkmark on the SpectraMax device icon. If instead there is a red circle with a slash through it, do the following:
 - a. Double-click the SpectraMax icon.
 - b. After the window pops up and populates with a list of instruments, click on the SpectraMax iD3 and connect to it. It should be the only instrument listed.
 - c. Close this pop-up window if it doesn’t automatically do so.
- 2B.5 Click on the orange gear icon to go to the Acquisition Setting menu.
- 2B.6 Choose/set the following:
 - a. Choose Fluorescence (not Abs or Lum).
 - b. Wavelengths: set to 480 and 530.
 - c. Plate type: 96 well, opaque.
 - d. Select all wells to read the entire plate, or you can omit empty cells/columns.
 - e. You can ignore all other settings (leave them as defaults).
- 2B.7 Open the plate reader drawer by clicking on the “Open/Close” button in the menu bar of the software, or by clicking the eject button on the iD3 touchscreen.
- 2B.8 Insert your plate, with A1 on the upper left and H12 on the bottom right.



- 2B.9 Close the drawer by clicking on the “Open/Close” button in the menu bar of the software, or by clicking the eject button on the iD3 touchscreen.
- a. DO NOT MANUALLY PUSH IN THE DRAWER BY HAND!!!
- 2B.10 Click the “Read” button (green arrow).
- 2B.11 Click “skip/okay” on the popup about optimization (you don’t need to check any boxes).
- 2B.12 The read will take ~1 minute 30 seconds.
- 2B.13 Once the run finishes, click “save” and save your experiment in your documents folder as a .sda file.
- a. This will allow you to reopen your data in SoftMax Pro if needed.
- 2B.14 Export your data in Excel (.xls) format:
- a. Click on the Plate icon in the top left of the screen.
 - b. Select Export.
 - c. Select your entire experiment (or plates of interest) & export your raw data.
 - d. Select “column” or “plate” format, depending on your preference.
 - e. Change file type to .xls.
- 2B.15 Open your file in Excel.
- 2B.16 Create a standard curve by creating a scatter plot of your standard data. The x-axis should be total ng of DNA from your standards, and the y-axis is their fluorescent readings.
- a. Check your assay type to determine total ng. For the HS kit, your standards will range from 0 to 100 total ng. The broad range kit will range from 0 to 1000 ng. **Note, this is TOTAL NG, not ng/μL!**
- 2B.17 Add a linear trendline to your standard curve, and display the equation and R² value.
- 2B.18 Remove any erroneous standards from your standard curve.
- 2B.19 Set the y-intercept to be the fluorescent value of your 0 ng standard (standard 1).
- 2B.20 Use the equation to calculate the total ng in each of your sample wells. You will be solving for x.
- a. Your y-value is the fluorescent reading from the assay.
- 2B.21 Divide these values by the number of μL of DNA used for each sample to determine your original ng/μL concentration.
- a. For example, if you only used 1 μL of DNA for your samples, then the outputted value IS your ng/μL.
 - b. If you used 5 μL of DNA, then divide the values by 5 to get your original ng/μL concentration.